

BMPro™-W

Serum free medium for
Antibody/Recombinant protein production



Suspended
CHO Cells

- ☐ Hydrolysates free
- ☐ Serum component free
- ☐ Low protein contents

Serum free medium for Antibody/Recombinant protein production

BMPro™-W

BMPro-W is serum-free complete medium formulated for CHO cells thriving in suspension culture systems. Tailored for optimal antibody and recombinant protein production, this medium is devoid of serum, hydrolysates, and maintains a low protein content. To use, simply supplement with 8mM L-glutamine.



Proliferation of CHO cell in suspension culture

Fig.1 Cell proliferation curve





Fig.2 Viability



Code (CSTI / NIPRO)	Product Name	Description	Volume
2502P10 / 87-015	BMPro-W	Serum free media for antibody & recombinant protein production (Optimized for CHO cell in suspension culture system) L-glutamine free	1000mL (P)
Supplement			
Code (CSTI / NIPRO)	Product Name	Description	Volume
1104 / 87-885	200mM L-Glutamine solution	200 mM L-glutamine Solution	100mL (P)

(P) PET bottle

※This product is for research use only. It cannot be used for therapeutic or diagnostic purpose.



1-16-16 Seikaen, Aoba-ku, Sendai,
989-3162, Japan
URL: <http://www.cstimedia.com>
E-mail: csti-info@ml.nipro.co.jp



NIPRO CORPORATION
3-26, Senriokashinmachi, Settsu,
Osaka, 566-8510, Japan