

Artificial serum

Serum substitute

(Supplement for serum-free culture media)



Culture media
supplement

Artificial serum
Animal-free




- For human Mesenchymal Stem Cells, human Fibroblasts, human Lymphocytes and other cell lines.
- A small amount of Artificial serum, 5% or less, supports cell growth as same as FBS.
- No lot-to-lot variation

Artificial serum (Animal-free)



- Artificial serum is developed as a substitute of serum (e.g. FBS).
- Additional "Artificial serum" to serum-free media enhance cell growth.
- Recommended concentration is 5% or less. Optimal concentration should be determined for specific use and cell type.
- The components of Artificial serum are animal-free.



1. Adhesion cells (Fibroblast)



2. Suspension cells (DAUDI)



※ Culture medium: ALyS505N-0

Code (CSTI / NIPRO)	Product name	Description	Storage	Volume
A2G20P2CC / 87-082	Artificial serum (Animal-free)	Animal-free Supplement	below -20°C	20mL (P)

(P) PET bottle

※ For Research use only.